

MEDICAL GAS EQUIPMENT



INDEX:

1. **Heyer Outlets and Accessories**
2. **Heyer Outlet Duplicators**
3. **Oxygen Flow Meters**
4. **Medical Gas Pressure Test Gauges**
5. **Double Gauge Medical Gas Regulators**
6. **Compact Medical Gas Regulators**
7. **Parts and Accessories**
8. **Anti-Static Medical Tubing**
9. **Oetiker Clamps**
10. **Heyer Probes**
11. **Medical Gas Fittings and Couplings**



National Sales and Training

Contact: Michele Hiscox

Cell: +27 62 344 1987

E-Mail: michele@trinitymedical.co.za

Office Sales Support

Contact: Dee Andrews

Cell: +27 72 199 0307

E-Mail: sales@trinitymedical.co.za

Wellness Equipment Sales and
Training

Contact: Sevy Andrews

Cell: +27 72 364 3027

E-Mail: sevy@trinitymedical.co.za

Office Admin

Contact: Estie Coetzer

E-Mail:

admin@trinitymedical.co.za



1a. SA Heyer Wall Outlet Sockets – Female

Order No	Description	Image	Order No	Description
MT0071TMS	Oxygen (O ²) Front Body Only	 <p>Front Body Only</p>	MT0130TMS	Oxygen (O ²) Front & Back Body
MT0072TMS	Nitrous Oxide (N ² O) Front Body Only		MT0131TMS	Nitrous Oxide (N ² O) Front & Back Body
MT0073TMS	Medical Air Front Body Only	 <p>Front & Back Body</p>	MT0132TMS	Medical Air Front & Back Body
MT0074TMS	High Pressure Air Front Body Only		MT0133TMS	High Pressure Air Front & Back Body
MT0075TMS	Vacuum Front Body Only		MT0134TMS	Vacuum Front & Back Body
MT0125TMS	Entonox Front Body Only		MT0135TMS	Scavenging Front & Back Body
	www.trinitymedical.co.za		MT0136TMS	Carbon Dioxide (CO ²) Front & Back Body

1b. SA Heyer Complete Wall Box / Flush Mount

Order No	Description	Image	Order No	Description
MT0087TMS	Oxygen (O ²) Wall Box	 <p>Wall Box</p>	MT0266TMS	Oxygen (O ²) Flush Mount
MT0088TMS	Nitrous Oxide (N ² O) Wall Box		MT0267TMS	Nitrous Oxide (N ² O) Flush Mount
MT0089TMS	Medical Air Wall Box	 <p>Flush Mount</p>	MT0268TMS	Medical Air Flush Mount
MT0090TMS	High Pressure Air Wall Box		MT0269TMS	High Pressure Air Flush Mount
MT0091TMS	Vacuum Wall Box		MT0270TMS	Vacuum Flush Mount
MT0092TMS	Scavenging Wall Box		MT0271TMS	Scavenging Flush Mount
MT0096TMS	Carbon Dioxide (CO ²) Wall Box		MT0272TMS	Carbon Dioxide (CO ²) Flush Mount

1c. SA Heyer Covers for Wall Box / Flush Mount

Order No	Description	Image	Order No	Description
MT0140TMS	Oxygen (O ²) Wall Box Surface Mount Cover	 <p>Wall Box Surface Mount Cover</p>	MT0160TMS	Oxygen (O ²) Flush Mount Cover
MT0141TMS	Nitrous Oxide (N ² O) Wall Box Surface Mount Cover		MT0161TMS	Nitrous Oxide (N ² O) Flush Mount Cover
MT0142TMS	Medical Air Wall Box Surface Mount Cover	 <p>Flush Mount Cover</p>	MT0162TMS	Medical Air Flush Mount Cover
MT0143TMS	High Pressure Air Wall Box Surface Mount Cover		MT0163TMS	High Pressure Air Flush Mount Cover
MT0144TMS	Vacuum Wall Box Surface Mount Cover		MT0164TMS	Vacuum Flush Mount Cover
MT0145TMS	Scavenging Wall Box Surface Mount Cover		MT0165TMS	Scavenging Flush Mount Cover
MT0146TMS	Carbon Dioxide (CO ²) Wall Box Surface Mount Cover		MT0166TMS	Carbon Dioxide (CO ²) Flush Mount Cover

1d. Medical Gas Outlet Conversion Kits Dräger/ MECO/ Gabler to SA Standard



Order No	Description	Image
MT01438TMS	Oxygen (O ₂)	
MT0439TMS	Nitrous Oxide (N ₂ O)	
MT0440TMS	Medical Air	
MT0441TMS	High Pressure Air	
MT0442TMS	Vacuum	
MT0443TMS	Scavenging	
MT0444TMS	Carbon Dioxide (CO ₂)	
MT0445TMS	Entonox	

* Kit includes Adaptor Plate, New Front Body, Screws and O-rings.

1e. Heyer Outlet Accessories

Order No	Description	Image	Order No	Description	Image
MT0094TMS	Eminence Medical Back Plate		MT0101TMS	Eminence Female Heyer Service Kit	
MT0095TMS	Eminence Medical Cover Mounting Bracket		MT0127TMS	Eminence Medical Female Socket Main Spring	
MT0098TMS	Eminence Medical Flush Mount Galvanized Boxes				

2. Heyer Outlet Duplicators

Description	Stock # Y-Shape	Stock # T-Shape	Image Y-Shape	Image T-Shape
Oxygen (O ₂)	MT0170TMS	MT0396TMS		
Medical Air	MT0171TMS	MT0397TMS		
Nitrous Oxide (N ₂ O)	MT0172TMS	MT0398TMS		
Vacuum	MT0173TMS	MT0399TMS		
High Pressure Medical Air	MT0174TMS	MT0400TMS		
Entonox	MT0175TMS	MT0401TMS		

3a. Oxygen Flow Meters with SA Heyer

Description	Order #	Image	Order #	Image
Oxygen (O ²) Flow Meter "Click Type" Output Setting: 0.5, 1,2, 3, 4, 5, 6 ,8 10, 15, 25 /min with Heyer connection (Aluminium Body)	X0035TMS (Single)		X0036TMS (Double)	
Oxygen (O ²) Flow Meter with Heyer connection, Ball and Tube type (0-15 /min) Brass Body (Spare parts available)	X0038TMS (Single)		X0039TMS (Double)	
Extended Oxygen (O ²) Flow Meter with Heyer connection, Ball and Tube type (0-15 /min) Brass Body (Extended to pass railing)	X0040TMS (Single)		X0041TMS (Double)	
Oxygen (O ²) Flow Meter - Gabler rail mount, ball and tube type with Heyer connection (0-15 /min)	X0042TMS (Single)		X0043TMS (Double)	
Oxygen (O ²) Flow Meter – Quad Ball & Tube-Type (0-15 /min) with Gabler rail mount.	X0397TMS			



The images shown are for illustration purposes only


3b. Oxygen Flow Meters with SA Heyer Connection

Description	Order #	Image
Blender Oxygen (O ²) Flow Meter Ball and Tube (0-15 /min)	X0389TMS	
Oxygen (O ²) Flow meter "Click Type" with female Heyer combination T-Bar and Male Heyer connection Output Setting: 0.5, 1,2, 3, 4, 5, 6 ,8 10, 15, 25 /min (Aluminium Body)	X0414TTMS	
Oxygen (O ²) Flow meter Ball and tube (0-15 /min) with female Heyer combination T-Bar and Male Heyer connection	X0415TMS	


3C. Oxygen Flow Meters with SA Heyer Connection

Description	Order #	Image	Order #	Image
Oxygen (O ²) Flow Meter Neonatal "Click Type" with SA Heyer connection 0-5 /min) (Aluminium Body)	X0057TMS (Single)		X0058TMS (Double)	
Oxygen (O ²) Flow Meter "Click Adjust" Output: 0.5, 1,2, 3, 4, 5, 6 ,8 10, 15, 25 /min with 1/8 NPT fitting	X0060TMS (Single)		X0061TMS (Double)	
Oxygen (O ²) Neonatal Flow Meter - Gabler rail mount, Click type with SA Heyer connection (0-5 /min)	X0046TMS (Single)			


4. Medical Gas Pressure Test Gauges Ø60 mm

Stock#	Gas Type	Image
MT0420TMS	Oxygen (O ²)	
MT0421TMS	Medical Air	
MT0422TMS	Vacuum	
MT0423TMS	High Pressure Air	
MT0424TMS	Nitrous Oxide (N ² O)	
MT0425TMS	Entonox	
MT0426TMS	Carbon Dioxide (CO ²)	

5a. Double Gauge Medical Gas Regulators

Description	Adjustable Flow Rate 0-15 /min		Adjustable Pressure 0-600 kPa		Image
	Pin Index Stock #	Bull Nose Stock #	Pin Index Stock #	Bull Nose Stock #	
Oxygen (O ²)	X0400TMS	X0406TMS	X0244TMS	X0382TMS	
Medical Air	X0401TMS	X0407TMS	X0245TMS	X0383TMS	
Nitrous Oxide (N ² O)	X0402TMS	X0408TMS	X0246TMS	X0384TMS	
Carbon Dioxide (CO ²)	X0403TMS	X0409TMS	X0247TMS	X0385TMS	
Entonox	X0404TMS	X0410TMS	X0248TMS	X0486TMS	





5b. Double Gauge Medical Gas Regulators with female Heyer

Description	Adjustable Pressure 0-600 kPa		Image
	Pin Index Stock #	Bull Nose Stock #	
Oxygen (O ²)	X0418TMS	X0428TMS	
Medical Air	X0419TMS	X0429TMS	
Nitrous Oxide (N ² O)	X0420TMS	X0430TMS	
Carbon Dioxide (CO ²)	X0421TMS	X0431TMS	
Entonox	X0422TMS	X0432TMS	








The images shown are for illustration purposes only






6a. Compact Medical Gas Regulators – Oxygen

Stock#	Description	Image
X0150TMS	Oxygen Regulator Pin-Index Economy version with brass chamber Input Pressure: 15100 kPa / 151 Bar / 2190 psi Output: 0, 0.5, 1, 2, 4, 5, 6, 8, 10, 15, 20 and 25 /min Click adjust with Barb connection 1 Year Warranty	
X0167TMS	Oxygen Regulator Pin-Index Superior Quality Input Pressure: 30000 kPa / 300 Bar / 4351 psi Output: 0, 0.5, 1, 2, 3, 4, 5, 6, 8, 10, 15, 20, 25 and 30 /min, Click Adjustment with DISS connection 2 Year Warranty	
X0151TMS	Oxygen Regulator Bull Nose Economy version with brass chamber Input Pressure: 15100 kPa / 151 Bar / 2190 psi Output: 0, 0.5, 1, 2, 4, 5, 6, 8, 10, 15, 20 and 25 /min Click adjust with Barb connection 1 Year Warranty	
X0168TMS	Oxygen Regulator Bull Nose Superior Quality Input Pressure: 30000 kPa / 300 Bar / 4351 psi Output: 0, 0.5, 1, 2, 3, 4, 5, 6, 8, 10, 15, 20, 25 and 30 /min, Click Adjustment with DISS connection 2 Year Warranty	





6b. Compact Medical Gas Regulators – Oxygen

Stock#	Description	Image
X0160TMS	Oxygen Regulator Pin-Index Economy version with brass chamber Input Pressure: 15100 kPa / 151 Bar / 2190 psi Output: 0, 0.5, 1, 2, 4, 5, 6, 8, 10, 15, 20 and 25 /min Click adjust with DISS connection 1 Year Warranty	
X0161TMS	Oxygen Regulator Bull Nose Economy version with brass chamber Input Pressure: 15100 kPa / 151 Bar / 2190 psi Output: 0, 0.5, 1, 2, 4, 5, 6, 8, 10, 15, 20 and 25 /min Click adjust with DISS connection 1 Year Warranty	
X0169TMS	Oxygen Regulator Pin-Index Superior Quality Input Pressure: 30000 kPa / 300 Bar / 4351 psi Output: 400 kPa Constant Pressure with DISS connection 2 Year Warranty	
X0170TMS	Oxygen Regulator Bull Nose Superior Quality Input Pressure: 30000 kPa / 300 Bar / 4351 psi Output: 400 kPa Constant Pressure, DISS connection 2 Year Warranty	
X0241TMS	Oxygen Regulator Pin-Index Superior Quality Input Pressure: 30000 kPa / 300 Bar / 4351 psi Dual Output: a) 400 kPa Constant Pressure, female SA Heyer connection b) 0,0.5, 1, 2, 3, 4, 5, 6, 8, 10, 15, 20, 25 and 30 /min, Click Adjustment, DISS connection 2 Year Warranty	

6c. Compact Medical Gas Regulators – Oxygen

Stock#	Description	Image
X0242TMS	Oxygen Regulator Bull Nose Superior Quality Input Pressure: 30000 kPa / 300 Bar / 4351 psi Dual Output: a) 400 kPa Constant Pressure, female SA Heyer connection b) 0,0.5, 1, 2, 3, 4, 5, 6, 8, 10, 15, 20, 25 and 30 /min, Click Adjustment, DISS connection 2 Year Warranty	
X0182TMS	Oxygen Regulator Pin-Index Economy version with brass chamber Input Pressure: 15100 kPa / 151 Bar / 2190 psi Dual Output: a) 0, 0.5, 1, 2, 4, 5, 6, 8, 10, 15, 20 and 25 /min Click adjust with Barb connection b) 414 kPa, constant pressure with female SA Heyer connection 1 Year Warranty	
X0183TMS	Oxygen Regulator Bull-Nose Economy version with brass chamber Input Pressure: 15100 kPa / 151 Bar / 2190 psi Dual Output: a) 0, 0.5, 1, 2, 4, 5, 6, 8, 10, 15, 20 and 25 /min Click adjust with Barb connection b) 414 kPa, constant pressure with female SA Heyer connection 1 Year Warranty	
X0174TMS	Oxygen Regulator Pin-Index Economy version with brass chamber Input Pressure: 15100 kPa / 151 Bar / 2190 psi Output: 0, 0.5, 1, 2, 4, 5, 6, 8, 10, 15, 20 and 25 /min Click adjust with DISS connection 1 Year Warranty	
X0175TMS	Oxygen Regulator Bull Nose Economy version with brass chamber Input Pressure: 15100 kPa / 151 Bar / 2190 psi Output: 0, 0.5, 1, 2, 4, 5, 6, 8, 10, 15, 20 and 25 /min Click adjust with DISS connection 1 Year Warranty	


6d. Compact Medical Gas Regulators – Oxygen

Stock#	Description	Image
X0184TMS	Oxygen Regulator Pin-Index Economy version with brass chamber Input Pressure: 15100 kPa / 151 Bar / 2190 psi Dual Output: a) 0, 0.5, 1, 2, 4, 5, 6, 8, 10, 15, 20 and 25 /min Click adjust with Barb connection b) 414 kPa, constant pressure with DISS connection 1 Year Warranty	
X0185TMS	Oxygen Regulator Bull-Nose Economy version with brass chamber Input Pressure: 15100 kPa / 151 Bar / 2190 psi Dual Output: a) 0, 0.5, 1, 2, 4, 5, 6, 8, 10, 15, 20 and 25 /min Click adjust with Barb connection b) 414 kPa, constant pressure with DISS connection 1 Year Warranty	
X0171TMS	Oxygen Regulator Bull Nose Superior Quality Input Pressure: 30000 kPa / 300 Bar / 4351 psi Dual Output: a) 400 kPa Constant Pressure, DISS Connection b) 0, 0.5, 1, 2, 3, 4, 5, 6, 8, 10, 15, 20, 25 and 30 /min, Click Adjustment, DISS connection 2 Year Warranty	
X0172TMS	Oxygen Regulator Pin-Index Superior Quality Input Pressure: 30000 kPa / 300 Bar / 4351 psi Dual Output: a) 400 kPa Constant Pressure, DISS connection b) 0,0.5, 1, 2, 3, 4, 5, 6, 8, 10, 15, 20, 25 and 30 /min, Click Adjustment, DISS connection 1 Year Warranty	

6e. Compact Medical Gas Regulators - Medical Air

Stock #	Description	Image
X0196TMS	Medical Air Regulator Pin Index Economy version with brass chamber Input Pressure: 15100 kPa / 151 Bar / 2190 psi Dual Output: a) 0, 0.5, 1, 2, 4, 5, 6, 8, 10, 15, 20 and 25 /min Click adjust with Barb connection b) 414 kPa, constant pressure with screw connection 1 Year Warranty	
X0197TMS	Medical Air Regulator Pin-Index Economy version with brass chamber Input Pressure: 15100 kPa / 151 Bar / 2190 psi Dual Output: a) 0, 0.5, 1, 2, 4, 5, 6, 8, 10, 15, 20 and 25 /min Click adjust with Barb connection b) 414 kPa, constant pressure with female SA Heyer 1 Year Warranty	
X0198TMS	Medical Air Regulator Bull-Nose Economy version with brass chamber Input Pressure: 15100 kPa / 151 Bar / 2190 psi Dual Output: a) 0, 0.5, 1, 2, 4, 5, 6, 8, 10, 15, 20 and 25 /min Click adjust with Barb connection b) 414 kPa, constant pressure with female SA Heyer 1 Year Warranty	
X0199TMS	Medical Air Regulator Bull-Nose Economy version with brass chamber Input Pressure: 15100 kPa / 151 Bar / 2190 psi Dual Output: a) 0, 0.5, 1, 2, 4, 5, 6, 8, 10, 15, 20 and 25 /min Click adjust with Barb connection b) 414 kPa, constant pressure with screw connection 1 Year Warranty	

6f. Compact Medical Gas Regulators – Other Gasses

Stock #	Description	Image
X0380TMS	CO2 Regulator Pin Index Economy version with brass chamber Input Pressure: 15100 kPa / 151 Bar / 2190 psi Output: 414 kPa, constant pressure with female SA Heyer 1 Year Warranty	
X0381TMS	N2O Regulator Pin Index Economy version with brass chamber Input Pressure: 15100 kPa / 151 Bar / 2190 psi Output: 414 kPa, constant pressure with female SA Heyer 1 Year Warranty	



The images shown are for illustration purposes only

7. Parts and Accessories


Order #	Description	Image	Order #	Description	Image
X0299TMS	O2 Leak Detection Concentrate 5 Litre		X0329TMS	O2 Leak Detection Spray Bottle	
X0354TMS	O2 Leak Detection Concentrate 1 Litre		X0053TMS	O2 Flow meter Control Knob	
X0052TMS	O2 Flow meter Spare Part Kit: Outer dome, Inner tube, Flow ball, O-Rings		X0281TMS	Adult Ventilator Test Lung 22mm Connection Autoclavable at 138°C	
X0054TMS	Humidifier bottle Reusable for O2 flow meter		X0283TMS	Neonatal Ventilator Test Lung Silicone Material 15mm Connection Autoclavable at	
X0178TMS	Oxygen female SA Heyer connection with 1/8 NPT thread		X0179TMS	Medical Air (LP) female SA Heyer connection with 1/8 NPT thread	
X0022TMS	X-Mas Tree Connector for Oxygen flow meter (Tail piece) Swivel type		X0051TMS	Metal Tail Piece Connector for Oxygen flow meter	
X0388TMS	Bodok (Pin-Index) Seals for O2, N2O and CO2		X0166TMS	O-Rings for Bullnose regulators	
MT0450TMS	Female Heyer Remote Extension on Gabler rail mount. ***		X0396TMS	Oxygen Regulator Female Heyer Extension Kit ***	
X0251TMS	Oxygen Analyser				

***Female Heyer remote extension and Female Heyer extension kit available in other gasses on request.




The images shown are for illustration purposes only



8a. Medical Tubing - Anti-Static Hoses

Order #	Description	Image	Order #	Description
MT0001TMS	Blue Anti-Static Medical Hose flexible 6 mm I.D (N2O) Pre-cut (Sold per 5 meter)		MT0020TMS	Blue Anti-Static Medical Hose flexible 6 mm I.D (N2O) 100-meter Coil
MT0002TMS	Black Anti-Static Medical Hose flexible 6 mm I.D. (Medical Air) Pre-cut (Sold per 5 meter)		MT0021TMS	Black Anti-Static Medical Hose flexible 6 mm I.D. (Medical Air) 100-meter Coil
MT0003TMS	White Anti-Static Medical Hose flexible 6 mm I.D. (Oxygen) Pre-cut (Sold per 5 meter)		MT0022TMS	White Anti-Static Medical Hose flexible 6 mm I.D. (Oxygen) 100-meter Coil
MT0004TMS	Yellow Anti-Static Medical Hose flexible 6 mm I.D. (Vacuum) Pre-cut (Sold per 5 meter)		MT0023TMS	Yellow Anti-Static Medical Hose flexible 6 mm I.D. (Vacuum) 100-meter Coil
MT0005TMS	Green Anti-Static Medical Hose flexible 6 mm I.D. (CO2) Pre-cut (Price per 5 meter)			



8b. Medical Tubing – Other

Order #	Description	Image	Order #	Description
MT0051TMS	Silicone Tubing Id: 6 mm and od 11 mm (1 Unit = Roll of 15 meters)		MT0050TMS	Silicone Tubing ID: 8 mm and OD-12 mm (1 Unit = Roll of 15 Meters)
MT0157TMS	Clear Tubing ID 8 mm (1 Unit = Roll of 30 Meter Roll)			

9. Oetiker Hose Clamps and Accessories

Order #	Description	Image	Order #	Description
MT0035TMS	OETIKER CLIP 2 Ear 7-9mm (1 Unit = 1 Packet of 10 Clamps)		MT0260TMS	OETIKER CLIP 2 Ear 7-9mm (1 Unit = 1 Packet of 100 Clamps)
MT0030TMS	OETIKER CLIP 2 Ear 11-13mm (1 Unit = 1 Packet of 10 Clamps)		MT0261TMS	OETIKER CLIP 2 Ear 11-13mm (1 Unit = 1 Packet of 100 Clamps)
MT0031TMS	OETIKER CLIP 2 Ear 13-15mm (1 Unit = 1 Packet of 10 Clamps)		MT0262TMS	OETIKER CLIP 2 Ear 13-15mm (1 Unit = 1 Packet of 100 Clamps)
MT0032TMS	OETIKER CLIP 2 Ear 15-18mm (1 Unit = 1 Packet of 10 Clamps)		MT0263TMS	OETIKER CLIP 2 Ear 15-18mm (1 Unit = 1 Packet of 100 Clamps)
MT0034TMS	OETIKER CLIP 2 Ear 17-20mm (1 Unit = 1 Packet of 10 Clamps)		MT0264TMS	OETIKER CLIP 2 Ear 17-20mm (1 Unit = 1 Packet of 100 Clamps)
MT0041TMS	OETIKER CLIP 2 Ear 20-23mm (1 Unit = Packet of 10 clamps)		MT0265TMS	OETIKER CLIP 2 Ear 20-23mm (1 Unit = 1 Packet of 100 Clamps)
MT0033TMS	OETIKER CLIP Pliers Side Pinch			









10. SA Heyer Probes – Male with Flow Restrictor Ball

Gas Type	6mm Barb	1/8 NPT thread	Image
Oxygen (O ₂)	MT0065TMS	MT0093TMS	 
Nitrous Oxide (N ₂ O)	MT0066TMS	MT0102TMS	
Medical Air	MT0067TMS	MT0103TMS	
High Pressure Medical Air	MT0068TMS	MT0104TMS	
Vacuum	MT0069TMS	MT0105TMS	
Entonox	MT0082TMS	MT0108TMS	
Carbon Dioxide (CO ₂)	MT0083TMS	MT0109TMS	

(Probes for other gasses available on request)
















www.trinitymedical.co.za

11. Medical Gas Fittings and Couplings

Stock #	Description	Image	Stock #	Description	Image
MT0190TMS	DISS O2 Female to 6mm Barb (O2 Blender fitting)		MT0191TMS	DISS Medical Air Male to 6mm Barb - Hex nut included	
MT0240TMS	DISS O2 Female to 6mm Barb		MT0230TMS	DISS O2 Hex Nut to 1/4" NPT Male	
MT0292TMS	DISS Male to 1/8 NPT Male with O-ring 1.5"		MT0234TMS	DISS O2 Body Adaptor to 1/4" NPT Male	
MT0291TMS	DISS Male to 1/8 NPT Male with O-ring		MT0235TMS	DISS Medical Air Body Adaptor to 1/4" NPT Male	
MT0228TMS	DISS O2 90 Degree Male to Female (Short)		MT0229TMS	DISS O2 90 Degree Male to Female (Long)	
MT0120TMS	DISS with nut to 6mm Barb to		MT0233TMS	DISS Cap & Chain Assembly	
MT0193TMS	NIST Medical Air Male to 6mm Barb		MT0237TMS	NIST Medical Air Female to 1/4" NPT Male	
MT0198TMS	NIST N2O Male to 6mm Barb		MT0239TMS	NIST N2O Female to 1/4" NPT Male	
MT0192TMS	NIST O2 Male to 6mm Barb		MT0236TMS	NIST O2 Female to 1/4" NPT Male	
MT0199TMS	NIST Vacuum Male to 6mm Barb		MT0238TMS	NIST Vacuum Female to 1/4" NPT Male	
MT0195TMS	NIST Medical Air Female to 6mm		MT0284TMS	Med Air Pin Index Yoke with	
MT0330TMS	NIST N2O Female to 6mm		MT0286TMS	N2O Pin Index Yoke with 1/4	
MT0194TMS	NIST O2 Female to 6mm Barb		MT0283TMS	O2 Pin Index Yoke with 1/4	
MT0329TMS	NIST VACUUM Female to 6mm		MT0285TMS	CO2 Pin Index Yoke with 1/4	

The images shown are for illustration purposes only

11. Medical Gas Fittings and Couplings

Stock #	Description	Image	Stock #	Description	Image
MT0287TMS	Stainless Steel Pigtail 36" (914mm) 1/4NPT to 1/4NPT 3000PSI		MT0293TMS	Elbow, 90Deg, 1/8 NPT Male x 1/8 NPT Female	
MT0288TMS	Stainless Steel Pigtail 48" (1219mm) 1/4NPT to 1/4NPT 3000PSI		MT0294TMS	Elbow, 90Deg, 1/8 NPT Female x 1/8 NPT Female	
MT0344TMS	1/4" NPT Male > 1/4 JIC Female Swivel (250Bar)		MT0227TMS	3-Way Connection	
MT0106TMS	Rectus 21 Connection Female 6mm Barb (Male = MT107)		MT0107TMS	Rectus 21 Connection Male 6mm Barb (Male = MT106)	
MT0124TMS	Scavenging Connection 19mm to 22mm barb with 16mm ID		MT0154TMS	Vacuum Heyer > 22mm Scavenge Fitting	
MT0123TMS	1/4" BSP Male to 6mm Barb		MT0128TMS	Brass T-Piece 6mm con	
MT0295TMS	Fitting 1/8 NPT Female to 6mm Hose Barb		MT0256TMS	Brass 1/8 NPT Plug	
MT0121TMS	6mm Barb to Male 1/8NPT thread		MT0448TMS	MC Suction Connection	



The images shown are for illustration purposes only