



ONEPRO Nebulizer

New Generation Mesh Mobile Nebulizer Nebuliser

- **High Therapeutic Quality:**

ONEPro is designed to ensure that you receive high therapeutic quality while going about your daily activities. This revolutionary device lets you **work, study, play, travel and move during your treatment.**

- **Silent and Portable:**

Combining the best features of ultrasound technology, ONEPro operates silently. Its compact and lightweight design offers maximum comfort during treatment.

- **Fast and Effective Nebulization:**

With compression technology, ONEPro guarantees quick and efficient nebulization of various medications (excluding hyaluronic acid). It is the ideal device for any type of use, offering versatility in managing your health.

- **Innovative “V” Design:**

ONEPro’s high-performance nebulization chamber features an innovative “V” design. This design ensures that medication is continually directed toward the center of the chamber, feeding the mesh continuously.

**EASIEST
NEBULIZER
EVER!**



LAPTOP POWERBANK SMARTPHONE

Technical Data:

- **Powered by:** USB port input 5V, 1A from external power feeding Cell phone, Tablet, Laptop, Power bank, Car USB socket.
- **Power Consumption:** 2.0 W
- **Neb-Rate:** 0.22 mL/min
- **MMAD (Mass Median Aerodynamic Diameter):** 1.51 μ
- **Respirable Fraction (< 5 μ m):** > 93%
- **Duty Cycle:** Non-stop
- **Noise Level:** \leq 35 dB (A)
- **Weight:** 18.0 g
- **Size (Main Unit):** 5.9 x 4.6 x 3.4 cm
- **Warranty:** 2 years